

Wheelchair cushions for medium-risk to high-risk patients



- Stable sitting ●
- Temperature reduction ●
- Prevents shearing ●
- Thermal dissipation ●
- Wound healing (SPPC and SPC) ●
- The only cushion on the market with a fan that supplies air ●



SPPC

High-risk Patients

The best positioning cushion with a fan.



SPC

High-risk Patients

The lower positioning cushion with fan, but with more pressure relief.



AireRx-cushions are built with an anatomically intelligent design to provide a stable sitting position with a superior pressure relief function. The cushions are designed to reduce heat, divert moisture that forms on the surface to provide maximum comfort and at the same time minimizing the risk of pressure sores. The cushion comes in three versions: SPPC, which is the most positioning cushion with a fan, SPC is less positioning but with more pressure relief and fan, GU is the easy and simple cushion without a fan which is ideal as basic cushion and for the active user.



AireRx
HEALTHCARE, LLC

AireRx's construction:


- Innovative anatomical design with leg separator.
- Built-in rechargeable battery-powered fan (for SPC & SPPC).
- Ventilation ducts for elimination of heat and moisture (for SPC & SPPC).
- Constructed of two types of pressure relieving foam.
- Maximum air permeable cover.
- The cushions are reusable and incontinence safe due to sealed surface.
- The cushion is equipped with carrying handles.



AireRx's properties:

- Very good ventilation thanks to the cushion's unique design and fan.
- The only cushion on the market with a fan that supplies air and ventilates!
- Very good wound-healing effect (for SPC & SPPC).
- Improved quality of life for the patient due to user can control the fan after needs. For example, affected moisture, heat and sweating.
- Better quality of life for all type of patients:
 - MS, which are usually warm.
 - With vaginal infections where fungus is counteracted due to the ability to control the temperature.
 - "wound users"
 - patients that are sitting for a long time.
- Clear labelling and easy positioning.
- Recommended maximum user weight approx 180 kg.
- 2 year warranty!

Sizes - see our website where they can be updated.
Washing instructions - see accompanying instructions.



Salubrious
Good Medical

